

die alternative*

THARMAC

CELLSTAIN - F

Automated Feulgenstainer for the DNA image cytometry

with

- 23 Position
- 1 slide holder for 30 or 50 slides
- 1 heating station, from 25°C to 28°C
- 1 fume extraction system

for 30 Slides: **Order-Nr. JC 600**

for 50 slides: **Order-Nr. JC 601**

Dimensions

L 160cm x H 45cm x D 55cm

Weight: approx. 38 kg

die alternative*

THARMAC

CELLSTAIN – F accessories optional

<u>Product</u>	<u>Order Nr.</u>
Feulgen Slide holder for 50 slides	JC 637
Feulgen Slide holder for 30 slides	JC 638
Slide dryer from 40°C to 60°C	JC 614
Running water station	JC 616
Programmer	JC 617
SmartCard	JC 618
trough 500 ml	JC 621
trough 750 ml (not for Xyleen)	JC 622
trough 750 ml (for Xyleen)	JC 623